

IN THE UNITED STATES PATENT AND TRADEMARK OFFICE

Art Unit : 1635  
Examiner : Louis V. Wollenberger  
Applicant : Nariyoshi Shinomiya et al.  
Appln. No. : 10/599,327  
Filing Date : March 12, 2007  
Confirmation No. : 7013  
For : *c-met* siRNA ADENOVIRUS VECTORS INHIBIT CANCER CELL  
GROWTH, INVASION AND TUMORIGENICITY

**DECLARATION UNDER 37 C.F.R. § 1.131**

We, the undersigned, do hereby declare as follows:

1. We are the co-inventors of the claims of the above-identified patent application.
2. The invention as defined in claims 1-12, 14-17, and 38 was conceived of and actually reduced to practice prior to January 6, 2003. The invention as defined in claims 13, 18-20, and 48-50 was conceived of and actually reduced to practice prior to July 7, 2003.
3. Evidence of our conception and reduction to practice of the invention as defined in claims 1-20, 38, and 48-50 is provided in the form of experimental data from the laboratory notebook of Nariyoshi Shinomiya, one of the named inventors (attached hereto as Exhibit A1-A10). More specifically, these laboratory notebooks show our development of an RNAi molecule directed to c-met:
  - a) in the cancer cell lines DU-145, SK-LMS-1, DA3, and M114 (Exhibit A1);
  - b) using siRNA expression vector pSilencer 1.0-U6 for human c-met (Exhibit A2 and A4);
  - c) targeting human c-met sequence 221, the target of SEQ ID No. 15 (Exhibit A3);
  - d) using pShuttle vector (Exhibit A5);

Applicant : Nariyoshi Shinomiya et al.  
For : *c-met* siRNA ADENOVIRUS VECTORS INHIBIT CANCER CELL  
GROWTH, INVASION AND TUMORIGENICITY  
Page : 2

e) using an Ad5 viral vector (Exhibit A6 and A9);

f) using a stable transformant (Exhibit A7 and A8); and

g) in DBTRG glioblastoma cells (Exhibit A10).

4. The documents attached as Exhibit A1-A10 were prepared contemporaneously with our conception and reduction to practice.


5. The acts referred to in the preceding paragraphs occurred in the United States.

6. The undersigned hereby declares that all statements made herein of his own knowledge are true and that all statements made on information and belief are believed to be true; and further, that these statements were made with the knowledge that willful false statements and the like so made are punishable by fine or imprisonment or both, under Sections 1001 of Title 18 of the United States Code and that such willful false statements may jeopardize the validity of the application or any patent issuing thereon.

September 30, 2009  
Date

September 26, 2009

Date

  
George F. Vande Woude

N. Shinomiya  
Nariyoshi Shinomiya

< Cell culture from the stock >

- ① DU-145 ——— RPMI 1640 + 5% FBS  
(-160°C stock = 1/11/02) (+ 1% pen/strep)
- ② SK-LMS-1 ——— DMEM + 10% FBS  
(passage 2) Arcc. (#11965-092) (+ 1% pen/strep)
- ③ DA3 ——— DMEM + 10% Calf Serum  
(-160°C stock = 3/23/01) Fetal Bovine Serum  
(+ 1% pen/strep)
- ④ M114 ——— DMEM + 10% Calf Serum  
(-160°C stock = 10-26-01) or Fetal Bovine Serum  
(+ 1% pen/strep)

Culture media & Cell preservation media

- \* 5% FBS - RPMI 1640 (Cat No: 11875-093  
Lot No: 1140467)
  - } FBS 15ml
  - } P/S 3ml
  - } RPMI 282ml

- \* 10% FBS - DMEM (Cat No: 11965-092  
Lot No: 1142807)
  - } FBS 30ml
  - } P/S 3ml
  - } DMEM 267ml

\* Cell preservation media

- |                                                                                                       |                                                                                                       |
|-------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|
| <ul style="list-style-type: none"> <li>} DMSO 7ml</li> <li>} FBS 20ml</li> <li>} RPMI 73ml</li> </ul> | <ul style="list-style-type: none"> <li>} DMSO 7ml</li> <li>} FBS 20ml</li> <li>} DMEM 73ml</li> </ul> |
|-------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------|

< Transformation of pSilencer<sup>TM</sup> 1.0-U6 vector  
into Top10 competent E. coli >

Top10 → from DeepFreezer

50  $\mu$ l → effendorf tube on ice

pSilencer 1.0-U6 0.2  $\mu$ l (= 150 ng)

↓  
25 min on ice

↓  
30 sec at 42°C

↓  
add 250  $\mu$ l SOC

↓  
shaking incubation (37°C, 250 rpm)

↓  
plate on amp LB plate (50  $\mu$ l, 100  $\mu$ l)



pSilencer 1.0-U6

< continued >

100  $\mu$ l → colonies are too many (can not be picked up)

50  $\mu$ l → ~ 2000 colonies

↓  
5 colonies are picked up



duplicated

plate ① → plasmid purification

② → stock of cells

↓  
37°C

# < Transformation mouse siRNA - part 4 >

## ① LB amp Agar

		pick up	
#57	8 colonies	(1, 2, 3)	
#60	1 colony	(1)	
#110	1 "	(1)	
#120	0 "	X	
#198	5 colonies	(1, 2, 3)	(1 = 0.5 hr)

LB amp medium 2ml  
shaking incubation

\* Transformation efficiency was not very good!

Since human siRNA oligo transformation obtained a good result, maybe there is a problem in the conc. of siRNA oligos.

Try again by using higher conc of siRNA oligos.

↳ #60, #110, #120 only

## ② Ligation again

multi pre Bbs I / Xba I digest & gel purified 0.5 µl

Xba I oligo 2 µl

5x Ligase buffer 2 µl

water 5 µl

T4 DNA Ligase 2.5 µl

Total 17 µl

Incubate at RT for 30 min

ADD EDTA  
50 µl 10 mM

Transformation  
in the same way as  
previously

4:25 PM

plasmid extraction using @Amax Mini-prep kit

(- 2.5 hr incubation)  
37°C, 250 strokes

Analysis by electrophoresis


(X) page 22b

# < Confirmation of the extracted plasmids >

## ① Digestion of plasmids with restriction enzymes

- #57, #178 → XbaI (Buffer 2)  
 #16-1, 2, 3, 4 #62-1, 2, 3, 4 #76-1, 2, 3, 4 #178-1, 2, 3, 4  
 → HindIII (Buffer 2)  
 \* 1ul each plasmid was digested at 37°C for 1hr

## ② Electrophoresis (1% agarose gel, GENE-MATE 150V for 1.5hr)



- 1 — λ HindIII
- 2 — mU6 (circular)
- 3 — #57
- 4 — #57/XbaI
- 5 — #178
- ⑥ — #178/XbaI
- 7 — λ HindIII
- 8 — pSilencer (circular)
- 9 — pSilencer/HindIII (2500)
- 10 — #16-1
- 11 — " 2
- ⑫ — " 3
- 13 — " 4
- 14 — #62-1
- ⑮ — " 2
- 16 — " 3
- 17 — " 4
- ⑱ — #76-1
- 19 — " 2
- 20 — " 3
- 21 — " 4
- 22 — #221-1
- ⑳ — " 2
- 24 — " 3
- 25 — " 4

\* gel stain  
 150ul water  
 + 30ul EtBr  
 (10ug/ml)

### (A) Mouse siRNA

XbaI digest → 3.4kb is OK

- (mU6-m443nt)  
 mU6 BbsI/XbaI ≈ 3.4kb  
 mU6 insert ≈ 3.4kb

### (B) human siRNA

after HindIII digestion  
 still circular is OK  
 (≈ 3.3Kb)

HindIII digest

○ → to be sequenced!

## ② Sequencing

#178 (#128-0) → mouse siRNA for e-Met (M13 R2 primer)  
#16 (#16-3) ]

16 (16-3)

#62 (#62-2)

第 24 号

#221 (#221)-2

→ human siRNA for c-Met (T3 primer (A))

Template plasmid: pUC19 (6.3 μg)

Primer 1.00 (2.00)

→ Sequence ord

~~7~~ 11/11/1913

M13 R2 = 66423 pool  $\rightarrow$  dissolved in 250  $\mu$ l of water  
(= 300 pool/ $\mu$ l)

T3 (A) : 2/140 pmol  $\rightarrow$  dissolved in 2 ml of water  $\rightarrow$  stock  
(200 pmol/ $\mu$ l) at -20°C

Then 10 fold dilution with water

priming water 10 gal + 10 gal + 250 gal used  
for equipping

Stock (250 pmol/ $\mu$ l) -20°C

### siRNA hairpin template sequences for human c-met RNAi

#### A. Criteria of Sequence Selection

1. 21mer that start with AA
2. GC content between 45-55%
3. No more than three consecutive T or G nucleotides can be present anywhere in the hairpin template sequences
4. The targeted region is selected from a given cDNA sequence beginning 50 to 100 nt downstream of the start codon. (5' or 3' UTRs and regions nearby the start codon are avoided, as these may be richer in regulatory protein binding sites. UTR-binding proteins and/or translation initiation complexes may interfere with binding of the siRNP or RISC endonuclease complex.)
5. Blasted selected sequences against Genbank
6. Nucleotide cross-match (no more than 17, much less would be better)
7. Start with AAG (for cloning into pSilencer 1.0-U6 vector system; RNA transcription begin with G, but according to the company's information a +1 G is not required.)
8. A pair of template sequences is as follows:

Sense Loop Antisense

Template 1: 5'-N(19)-TTCAAGAGA-N(19)-TTTTT-3' (53mer)

Template 2: 3'-CCGG-N(19)-AAGTTCTCT-N(19)-AAAAAATTAA-5' (61mer)

Apal

EcoRI

Target sequence (161): AAGACCTCAGAAGGTTGCTG

Blast search (In new window)

Position in gene sequence: 415

GC content: 47.6%

siRNA Sense strand: GACCTCAGAAGGTTGCTGTT

siRNA Antisense strand: CAGCAACCTTCTGAAGGTTCT

#16-1: 5'-GACCTCAGAAGGTTGCTG-TTCAAGAGA-CAGCAACCTTCTGAAGGTC-TTTTT-3'

#16-2: 3'-CCGG-CTGGAAGTCTTCCAACGAC-AAGTTCTCT-GTCGTTGGAAGACTTCCAG-AAAAAATTAA-5'

#16-2: 5'-AATTAATAAAA-GACCTCAGAAGGTTGCTG-TCTCTTGAA-CAGCAACCTTCTGAAGGTC-GGCC-3'

Target sequence (62): AAGCCAGATTCTGCCGAACCA

Blast search (In new window)

Position in gene sequence: 1236

GC content: 52.4%

siRNA Sense strand: GCCAGATTCTGCCGAACCAT

siRNA Antisense strand: TGGTTCGGCAGAACTCGGCTT

#62-1: 5'-GCCAGATTCTGCCGAACCA-TTCAAGAGA-TGGTTCGGCAGAACTCGGC-TTTTT-3'

#62-2: 3'-CCGG-CGGTCTAAGACGGCTTGGT-AAGTTCTCT-ACCAAGCCGCTTAGACCC-AAAAAATTAA-5'

#62-2: 5'-AATTAATAAAA-GCCAGATTCTGCCGAACCA-TCTCTTGAA-TGGTTCGGCAGAACTCGGC- GGCC-3'

Target sequence 76: AAGCGCGCGTGATGAATATC

Blast search (In new window)

Position in gene sequence: 1417

GC content: 52.4%

siRNA Sense strand: GCGCGCGTGATGAATATCTT

siRNA Antisense strand: GATATTCATCACGCGCGCTT

#76-1: 5'-GCGCGCGTGATGAATATC-TTCAAGAGA-GATATTCATCACGCGCGC-TTTTT-3'

#76-2: 3'-CCGG-CGCGCGGCACTACTATAG-AAGTTCTCT-CTATAAGTAGTGCCGCGCG-AAAAAATTAA-5'

#76-2: 5'-AATTAATAAAA-GCGCGCGTGATGAATATC-TCTCTTGAA-GATATTCATCACGCGCGC- GGCC-3'

Target sequence (221): AAGTCAGTATCCTCTGACAG

Blast search (In new window)

Position in gene sequence: 3310

GC content: 47.6%

siRNA Sense strand: GTGCAGTATCCTCTGACAGTT

siRNA Antisense strand: CTGTCAAGGATACCTGCACCT

#221-1: 5'-GTGCAGTATCCTCTGACAG-TTCAAGAGA-CTGTCAAGGATACCTGCAC-TTTTT-3'

#221-2: 3'-CCGG-CACGTCATAGGAGACTGTC-AAGTTCTCT-GACAGTCTCTATGACGTG-AAAAAATTAA-5'

#221-2: 5'-AATTAATAAAA-GTGCAGTATCCTCTGACAG-TCTCTTGAA-CTGTCAAGGATACCTGCAC- GGCC-3'



# < Sequence confirmation of siRNA plasmid >

## Expression plasmids for c-met RNAi (final clones)

### 1. Mouse siRNA (host plasmid = mU6pro)

Target sequence No.	Position in gene sequence	Ligation of synthesized oligos & transformation	Sequence confirmation of the inserted oligos	Final clone #s
#57	950	finished	confirmed	#57-1, #57-2, #57-3
#60	988	finished	confirmed	#60-1, #60-4
#110	1839	finished	confirmed	#110-1, #110-2, #110-3
#120	1977	Very low transformation efficiency	-	-
#178	2762	finished	confirmed	#178-0, #178-1, #178-2

### 2. Human siRNA (host plasmid = pSilencer)

Target sequence No.	Position in gene sequence	Ligation of synthesized oligos & transformation	Sequence confirmation of the inserted oligos	Plasmid amplification for transfection
#16	415	finished	confirmed	#16-3
#62	1236	finished	confirmed	#62-2
#76	1417	finished	Wrong sequence	-
#221	3310	finished	confirmed	#221-6

#red, #blue : plasmid clones were amplified and used for transfection

Control plasmid: mU6pro = mU6#1, pSilencer = pSil#1

< plasmid amplification > #110-1 #221-6 #mU6#1 pSil#1

2 ml → shaking incubation  
LB amp

#110-1  
#221-6 → from LB plate  
#mU6#1  
pSil#1 → from -80°C stock

100 ml overnight shaking incubation !!

EXHIBIT A4

## < Large Scale Plasmid Amplification

— Murine c-Het siRNA expression plasmids >

(1) Yesterday → separation culture from -80°C stock

(2) Today → pick up one colony & streaking culture

in 2 mL LB-Amp  
{ 100 µg/mL }  
[ 9:02 a.m. ]

### Ligation reaction (brief protocol)

#### a. Stick-ends (Cohesive Ends) Ligation

1. To an eppendorf tube add the following:

5x ligase reaction buffer	4 µL
vector DNA	3-30 fmol (2.5-25 ng)
insert DNA	9-90 fmol (7.5-75 ng)
(total DNA	0.01-0.1 µg)
autoclaved distilled water	to 19 µL
<u>T4 DNA Ligase</u>	<u>1 unit (in 1 µL)</u>

Final volume 20 µL

2. Mix gently. Centrifuge to bring the contents to the bottom of the tube.
3. Incubate at 23°C to 26°C for at least 5 min (30 min would be better).
4. Add 1 µL of 0.5 M EDTA to inactivate the enzyme.
5. Store the reaction at 4°C.
6. Dilute an aliquot of the ligation reaction five-fold in autoclaved distilled water and use it to transform competent cells.

#### b. Blunt ends Ligation

1. To an eppendorf tube add the following:

5x ligase reaction buffer	4 µL
vector DNA	15-60 fmol (25-250 ng)
insert DNA	45-180 fmol (75-750 ng)
(total DNA	0.1-1.0 µg)
autoclaved distilled water	to 19 µL
<u>T4 DNA Ligase</u>	<u>1 unit (in 1 µL)</u>

Final volume 20 µL

2. Mix gently. Centrifuge to bring the contents to the bottom of the tube.
3. Incubate at 14°C for 16-24 hr.
4. Add 1 µL of 0.5 M EDTA to inactivate the enzyme.
5. Store the reaction at 4°C.
6. Dilute an aliquot of the ligation reaction five-fold in autoclaved distilled water and use it to transform competent cells.

# Design (restriction enzyme digestion & ligation)

## a. Murine c-met siRNA

pShuttle (vector): XbaI/HindIII double digestion --> gel purification  
mU6pro (insert): HindIII/XbaI double digestion --> gel purification

XbaI = REact 2, 37°C, 1 hr  
HindIII = REact 2, 37°C, 1 hr


\*stick ends ligation (direction of the insert is reverse)

## b. Human c-met siRNA

pShuttle (vector): KpnI digestion --> phenol extraction/ethanol precipitation  
EcoRV digestion --> gel purification  
pSilencer (insert): KpnI/SmaI double digestion (30°C then 37°C) --> gel purification

KpnI = REact 4, 37°C, 1 hr  
SmaI = Gene Choice buffer 4 (=REact 4), 30°C, 1 hr  
EcoRV = Gene Choice buffer 2 (=REact 2), 37°C, 1 hr


\*blunt ends ligation (the same direction)



## &lt; Western blot ~ Results &gt;

① Control

② Hgfi-pretreated



\* 50-190 kDa → phosphorylation bands are too strong

↓  
 [ No time load 50µg protein  
 Separate with 8% gel

## &lt; Adenovirus purification using Virapur &gt;

	OD <sub>260</sub>	280	260/280
pAd ②	0.161	0.135	1.195
pAd ③	0.186	0.125	1.495

titre (cpu/ml)

 $1.8 \times 10^{12}$  $2.0 \times 10^{12}$ 

pAd ②  
 pAd ③ ] → harvest tomorrow

pAd ②, ③  
 300µl each × 9  
 + 1 each

< Western blot ~ siRNA transfectant

SK-LMS-1 & DU-145 >

primary Ab


( c-Met (C28) 1:2000  
β-Actin (AC15) 1:5000

protein

20 μg was loaded in each lane  
(both SK and DU)

secondary Ab

( α-rabbit 1:1000  
α-Mouse 1:5000



C-Met Downregulation is also observed in DU-145 cells

# < Histology of the remnant livers >

Serial number	Sample name	Serial number	Sample name
GW-4001	C 12h-1	GVW-4016	H 12h-1
4002	C 12h-2	4017	H 12h-2
4003	C 12h-3	4018	H 12h-3
4004	C 24h-1	4019	H 24h-1
4005	C 24h-2	4020	H 24h-2
4006	C 24h-3	4021	H 24h-3
4007	C 36h-1	4022	H 36h-1
4008	C 36h-2	4023	H 36h-2
4009	C 36h-3	4024	H 36h-3
4010	C 48h-1	4025	H 48h-1
4011	C 48h-2	4026	H 48h-2
4012	C 48h-3	4027	H 48h-3
4013	C 72h-1	4028	H 72h-1
4014	C 72h-2	4029	H 72h-2
4015	C 72h-3	4030	H 72h-3

Difference in the mitotic indices?



4001 ~ normal, hard to find mitotic cells

4002 ~ small vacuoles ⊕

4003 ~ bleeding inside the liver, hyaline-like changes

4004 ~ many vacuoles ⊕

4005 ~ small vacuoles ⊕  
(around the vein, no vacuoles)

4006 ~ "

4007 ~ vacuoles ⊕

4008 ~ small vacuoles ⊕

4009 ~ small vacuoles ⊕, normal

4010 ~ small vacuoles ⊕ / mitotic cells occasionally

4011 ~ vacuoles ⊕ / mitotic cells, occasionally

4012 ~ " / "

4013 ~ mitotic cells 1 ~ 3 / HPF

4014 ~ "

4015 ~ "

4016 ~ many bleeding sites

4017 ~ small vacuoles ⊕

4018 ~ " ⊕

4019 ~ vacuoles ⊕

4020 ~ " ⊕

4021 ~ " ⊕

4022 ~ vacuoles ⊕ depends on the portion  
Mitosis occasionally

4023 ~ vacuoles ⊕, Mitosis occasionally

4024 ~ vacuoles ⊕ hyaline necrosis  
inflammation

4025 ~ vacuoles ⊕ Mitosis 2 ~ 3 / HPF

4026 ~ mitosis 2 ~ 3 / HPF

4027 ~ "

4028 ~ Mitosis 2 ~ 3 / HPF

4029 ~ "

4030 ~ "

< Western blot ~ siRNA stable transformant >  
Expression of EGFR

□ SK-LMS-1

X, M, ①, 1-i, 2-i, 3-j, 4-e, 5-b, 5-c, 1-k, 4-f, X

□ DU-145


X, M, ①, 1-c, 1-d, 1-f, 2-l, 3-e, 3-h, 4-g, 4-h, X

200pg protein  
each  
loaded

1st Ab:  $\alpha$ -EGFR 1:1000

2nd Ab:  $\alpha$ -Rabbit 1:2000

Results ~



- ① { in SK-LMS-1 → expression of EGFR varies & also correlates well with proliferative activity  
in DU-145 → no remarkable changes in the EGFR levels

< Large scale virus production ~ pad ① >

T175 x 5 flasks



Virapour kit

OD<sub>260</sub>

OD<sub>280</sub>

0.087


0.056

$\frac{260}{280} = 1.554$

$$* \text{cpu/ml} = 0.087 \times 50 \times 1.1 \times 10^{12} = 4.8 \times 10^{12} (\text{cpu/ml})$$

500  $\mu$ l each x 8 tubes + 1 additional tube  $\rightarrow -80^{\circ}\text{C}$   
stored

< Cell scattering activity ~ B43 cells >



10 ng/ml HGF was added to each well  
( 9:30 a.m. ~ 3:00 p.m. )

< Large scale virus production ~ pAd ⑤ >

T175 x 5 flasks



Virapour kit

OD<sub>260</sub>

OD<sub>280</sub>

$\frac{260}{280}$  ratio

0.041

0.029

1.394

$$\text{particle number} = 0.041 \times 50 \times 1.1 \times 10^{12} = 2.3 \times 10^{12} (\text{cpu/ml})$$

500  $\mu$ l each x 8 tubes + 1 additional tube  $\rightarrow -80^{\circ}\text{C}$   
( > 1000  $\mu$ l ) stored



# < RT-PCR results >

sid cont      sid HQF-3

\* after 1/3 pH

0	1	2	3	4	5	6	7
0h	2h	6h	12h	24h	36h	48h	72h

← hHQF  
(539bp)


↑  
no hHQF expression in sid control mice.

← hHQF (365bp) → no remarkable changes in terms of time course.

← mβ-actin  
(376bp)

< DBTRG cells siRNA ads. infect  $\Rightarrow$  HGF stimulation >

## siRNA suppresses Met phosphorylation



- Met is down regulated in pAd9 infected cells
- Met phosphorylation is also significantly suppressed in pAd9 infected cells
- But regarding to the Akt-phosphorylation & p44/42 MAPK phosphorylation, strong phosphorylation bands were observed from the beginning (before HGF stimulation)